# ▶ THC BodyPlus

(D)	<i>Bedienungsanleitung – BodyPlus IR-Thermometer</i>	A - 1
(GB)	<i>Operating manual – BodyPlus IR-Thermometer</i>	B - 1
(F)	<i>Manuel d'utilisation – Thermomètre IR BodyPlus</i>	C - 1
(I)	<i>Istruzioni per l'uso – Termometro IR BodyPlus</i>	D - 1
(NL)	<i>Bedieningshandleiding – BodyPlus IR-Thermometer</i>	E - 1
(E)	<i>Manual de instrucciones – Termómetro IR BodyPlus</i>	F - 1
(P)	<i>Manual de instruções – Termómetro IR BodyPlus</i>	G - 1
(PL)	<i>Instrukcja obsługi – Termometr na podczerwień BodyPlus</i>	H - 1
(TR)	<i>Kullanım kılavuzu – BodyPlus IR Termometre</i>	I - 1
(RUS)	<i>Инструкция по эксплуатации – инфракрасный термометр</i>	J - 1
(DK)	<i>Bedjeningsvejledning – BodyPlus IR-termometer</i>	K - 1
(FIN)	<i>Käyttöohje – BodyPlus infrapunakuumemittari</i>	L - 1
(N)	<i>Bruksanvisning – BodyPlus IR-termometer</i>	M - 1
(S)	<i>Bruksanvisning – BodyPlus IR-termometer</i>	N - 1

CE 0197

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
This infrared thermometer has been built in accordance with the present state of the art. The device complies with the following standards: EN 60601-1-2:2007. Please contact our service technicians on +49 (0) 2452 962 530 if you have any questions regarding electromagnetic tolerances or the measured values. This conformity has been proven; the manufacturer is in possession of the relevant declarations and documents.

## 1. General description

This infrared thermometer for non-contact measurement has been specially developed in order to measure the temperature of persons irrespective of the ambient temperature. Temperature deviations may occur due to different types and thicknesses of skin.

In the case of a strongly differing room temperature, the thermometer should be placed in the room concerned 15 to 20 minutes before performing the first measurements.

Electronic equipment may not be disposed of as domestic waste, but must be correctly disposed of in accordance with the EUROPEAN PARLIAMENT AND EUROPEAN COUNCIL Directive 2002/96/EC of 27 January 2003 regarding used electrical and electronic equipment. Please dispose of this equipment in accordance with the applicable legal regulations when its service life is ended.

## 2. Safety instructions

- This appliance may only be used for the purposes described in the user manual.
- The appliance may only be used at ambient temperatures from 10 ° to 40 °C.
- The appliance should not be exposed to power surges.

- The thermometer should not be exposed to extreme temperatures below -20 °C or above +50 °C.
- The appliance may not be used if the relative humidity exceeds 85 %.
- The glass in front of the lens is cleaned with a cotton wool bud lightly moistened with a 70 % alcohol solution.
- If you should encounter problems with your new appliance, please contact the seller. Never attempt to repair it yourself.

### 3. Functions

- Precise non-contact measurements
- User-defined temperature selection (°C or °F)
- Freely selectable body and surface measurements
- Freely selectable alarm value
- Memory function for the last 32 measurements
- Automatic saving of the current measured value and automatic power-off
- Automatic selection range and display resolution 0.1 °C (0.1 °F)
- Backlit LCD display

### 4. Intended use

The non-contact infrared thermometer has been specially developed for non-contact temperature measurements on the surface of the body and on the forehead of small children and adults without touching the body of the person concerned. The thermometer is intended both for private use at home and as a reference appliance for doctors and clinics.

The appliance has been manufactured in compliance with the ASTM 1965-1998 standard specifications for infrared thermometers for the provisional determination of the temperature of patients.

In addition, the surface temperature function of the non-contact infrared thermometer can be used, for example, to determine the temperature of a baby's bottle or bath water. This conforms to the EMC directive 89/336/EEC.

The declaration of conformity in accordance with 92/43/EEC applies only to the medical use of the product. Use for all non-medical measurements - whether non-contact or by means of contact - was not a subject of the conformity appraisal procedure in accordance with the above-mentioned directive.

The appliance conforms to the specifications of EEC Directive 93/42 for medical products.

## 5. Precautions to be taken before use

- After moving the appliance to a new location with a considerable temperature difference to the previous location, it is strongly recommended to allow the non-contact infrared thermometer to acclimatise for 15 - 20 minutes in the new surroundings before taking the first measurements. Only in this way can you ensure that the same stable measurement results are obtained.
- Temperature deviations may occur due to different types and thicknesses of skin. The infrared thermometer should not be exposed to direct sunlight or water. Maintain a pause of approx. 1 minute between individual measurements. If you need to measure your temperature more frequently over a short period of time, the measurement results may fluctuate slightly.
- Do not measure your temperature immediately after eating, sports training or bathing. Wait at least half an hour.
- The temperature of the human body fluctuates within a natural range. There is no standard body temperature.
- Measure and note down your body temperature on normal days. This way you can recognise a raised temperature or fever more easily.

### Normal temperatures for the respective measurement methods:

Rektal	36,6 - 38,0 °C (97,8 - 100,4 °F)
Oral	35,5 - 37,5 °C (95,9 - 99,5 °F)
Armpit	34,7 - 37,3 °C (94,4 - 99,1 °F)
Ear	35,8 - 38,0 °C (96,4 - 100,4 °F)
Temple (Thermo Flash)	35,8 - 37,8 °C (96,4 - 100,0 °F)

A person's body temperature can vary over the course of the day. It can additionally be affected by other factors, such as age, sex, type of skin or thickness of skin.


### Normal temperatures based on age:

0 - 2 years	36,4 - 38,0 °C (97,5 - 100,4 °F)
3 - 10 years	36,1 - 37,8 °C (97,0 - 100,0 °F)
11 - 65 years	35,9 - 37,6 °C (96,6 - 99,7 °F)
> 65 years	35,8 - 37,5 °C (96,4 - 99,5 °F)

**⚠ The use of this non-contact infrared thermometer does not in any way replace medical treatment.**

## 6. Appliance illustration

1. IR-Sensor
2. LCD display
3. Mode selection
4. Button ▼
5. Button ▲
6. MODE button
7. Trigger (starts the measurement)
8. Battery compartment lid

## 7. Operation

1. The correct function of the non-contact infrared thermometer should be checked before using it for the first time by following the steps below. In order to ensure that the thermometer always delivers reliable and stable temperature measurement results, measure the temperature of a test person using a conventional thermometer. In a normal case, this may be, for example, 37.5 °C (99.5 °F).
2. Now measure the same test person's temperature with the non-contact infrared thermometer, maintaining a distance of between approx. 5 cm and maximum 15 cm between the thermometer and the forehead of the person concerned. When doing so, care should be taken that the measurement result is not distorted, for example, by hair hanging down or beads of sweat. Beads of sweat must be wiped away before measurement. If the value measured by the non-contact infrared thermometer is also approx. 37.5 °C (99.5 °F), i.e. it corresponds to the value obtained with the conventional thermometer, then the settings are correct and the infrared thermometer can be used immediately without problem.

If you obtain a lower measured value using the non-contact infrared thermometer, e.g. 36.4 °C (97.4 °F), then this value lies in the above example exactly 1.1 °C (2.2 °F) below the guide value of the conventional thermometer. The tempera-

ture measured by the non-contact IR thermometer must then be corrected upwards by 1.1 °C (2.2 °F). To do this, press and hold the MODE button for 2 seconds until the combination F1 appears on the display. Now press the same button again; F3 is displayed. Then press the ▲, button and enter the value of 1.1 °C (2.2 °F) given in the example.

3. In order to ensure that the steps listed above have been performed correctly and that you get the maximum possible benefit from your new non-contact infrared thermometer, we recommend that you take the test person's temperature once again with the infrared thermometer.

## 8. Taking measurements

First of all, insert the batteries in the appliance, ensuring the correct polarity.

After inserting the batteries, the appliance requires a warm-up period of around 10 minutes before it can be put into operation. This also applies to all further battery replacements.

Aim at the person's forehead, maintaining a distance of approx. 5 cm (2.5 inches). Now press the trigger. The temperature will be displayed immediately.

Before measuring, make sure that the forehead is not covered by hair hanging down, beads of sweat or other objects, which could influence the measurement. It is therefore essential to

wipe away beads of sweat before every measurement. The measurement can optionally also be taken behind the ear lobe.

### Temperature scale setting – F1

Press and hold the MODE button for approx. 2 seconds. F1 appears on the display. Press the ▼ button in order to display the temperature in Celsius and the ▲ button to display it in Fahrenheit.

### Alarm function – F2

This mode can be used to compensate for temperature differences in your non-contact infrared thermometer. Press and hold the MODE button for 2 seconds. F1 appears on the display. Press the MODE button twice to access F3. Press the ▲ to increase the limit value and ▼ to lower the threshold value. The set value is around 38.0 °C (100.4 °F).

### Temperature difference setting – F3

This mode can be used to compensate for temperature differences in your non-contact infrared thermometer. Press and hold the MODE button for 2 seconds. F1 appears on the display. Press the MODE button twice to access F3. Press the ▲ button to increase the temperature and the ▼ button to decrease it. It is essential to perform this check if temperature differences arise due to seasonal or changed ambient conditions. This applies only to Body mode.

## Alarm ON/OFF – F4

Press and hold the MODE button for approx. 2 seconds. F1 appears on the display. Press the MODE button three times to access F4. Press the ▲ button to activate the alarm. A short acoustic signal sounds and the  symbol appears on the display. Press the ▼ button to deactivate the alarm. 'OFF' appears on the display.

## Exiting mode setting

Press and hold the MODE button until the display switches off or wait a few seconds until the appliance switches itself off. The non-contact infrared thermometer has been specially developed to measure body temperature. Use the Body mode for this. The measurement range for this mode is from 32.0 °C to 42.0 °C (86.0 °F – 108.0 °F).

The non-contact infrared thermometer can also be used to measure the temperature of objects, e.g. (baby) food or liquids (baby's bottles), or other surfaces. Use the Surface mode for this. The measurement range for this mode is from 0.0 °C to 60.0 °C (32.0 °F – 140.0 °F).

 ***The temperature at the surface of the body always differs from the temperature inside the body. Therefore, always use Body mode to measure the internal body temperature.***

Also make sure that you have selected Surface mode if you wish to determine the surface temperature of an object.

## Saving of data

The temperature measurement is saved automatically. The value appears at the top right-hand corner on the LCD display. Press ▲ or ▼ to read off the last temperature measurement. The ▲ and ▼ buttons must be pressed simultaneously for 2 seconds with the appliance switched off in order to view the last temperature measurement. If you set the order to 0 in this setting and then press the MODE button, all previous temperature measurements will be deleted.

## Battery replacement

If the symbol  lights up on the display, this means that the batteries are spent. In this case, proceed as follows: Open the lid of the battery compartment and replace the spent batteries with new ones, ensuring that the polarity of the batteries is correct. Inserting the batteries incorrectly could damage the appliance and result in any warranty claims being voided. Do not use rechargeable batteries. Use only high-grade, non-rechargeable batteries.

## Lifetime

The non-contact infrared thermometer has been developed to meet the highest demands in daily and professional use. The lifetime is around 6 Years.

## 9. Cleaning and Maintenance

### Cleaning and disinfection

Please check the housing of your infrared thermometer at regular intervals for any visible signs of damage or soiling and please make sure that the keys are in good working order. Please ensure that the infrared lens of the infrared sensor is always clean and free of any dirt or other substances. Your infrared thermometer should be checked once a year using a suitable means of testing to check the measuring accuracy (see operating instructions). The monitoring of the measuring accuracy, which must take place at least once a year, is the responsibility of the user. The manufacturer, which has suitable normative values at its disposal, can be commissioned to perform comparative measurements.

Please use a slightly moistened soft cloth and a soap solution to remove any dirt or other signs of fouling on the housing or the infrared lens. Please make sure that there are no remains of the cleaning solution on either the housing or the lens when cleaning has been completed. Do not immerse the thermometer in any liquid! The thermometer can then be disinfected using a 70% isopropyl alcohol solution. Please store the thermometer at room temperature in its original casing or a protective sheath when not used over a longer period.

## 10. General information

- The most important and at the same time most fragile part of your infrared thermometer is the glass that protects the lens. Please be very careful with this.
- The glass is cleaned with a cotton cloth slightly moistened with water or a 70 % alcohol solution.
- Never use batteries other than the AA batteries described above. Do not use rechargeable batteries. Never throw spent batteries on an open fire.
- Remove the batteries from the infrared thermometer if the appliance is not to be used for a longer period of time.
- Do not expose the IR thermometer to sunlight for long periods, to direct sunlight, or to water.
- Protect the appliance against shocks. Even slight jolts can damage the measuring instrument.

## 11. Troubleshooting

The following points are intended to help you to solve certain problems that may occur when using the non-contact infrared thermometer.

If the problem cannot be solved by means of the following points, please contact Customer Service.

**The temperature shown on the display exceeds 60 °C.**

- The temperature is displayed in Fahrenheit. Change the temperature setting to Celsius by pressing the MODE button to access F1.

**The body temperature shown on the display is lower than 32 °C (89.6 °F).**

- The appliance is in Surface mode. The temperature shown on the display indicates the temperature radiated by the surface of the skin. Switch to Body mode if necessary.

**The message 'Hi' appears on the display.**

- The temperature analysis lies outside the selected measurement range, i.e. above 42.5 °C (108 °F) in Body mode or above 60 °C (140 °F) in Surface mode.

**The message 'Lo' appears on the display.**

- The temperature analysis lies outside the selected measurement range, i.e. below 32.0 °C (90 °F) in Body mode or below 0 °C (32 °F) in Surface mode.
- The measurement cannot be performed due to hair or beads of sweat on the skin. Make sure that no hair hanging down or beads of sweat are preventing correct measurement .
- An air current is preventing measurement. Make sure that no draughts can impair the sensitive infrared measurement.
- The time interval between two measurements is too small. Make sure that there is a pause of at least 15 seconds between two measurements – the recommended pause between two measurements is one minute.
- The measurement distance is too large. Please ensure that you

maintain the recommended measuring distance of between 5 and 15 cm (2.5 and 5.9 inches).

## 12. Technical data

### Operating conditions

Display resolution. . . . .	0,1 °C (0,1 °F)
Operating temperature. . . . .	10 bis 40 °C (50 bis 140 °F)
Humidity . . . . .	≤ 85%
Power supply . . . . .	DC 3V (2x AA batteries)
Dimensions . . . . .	149 x 77 x 43mm / 5,9 x 3 x 1,7 Zoll (LxBxH)
Weight . . . . .	Overall weight 400 g / net weight 172 g

### Measurement range

for body measurements. . . . .	32,0 to 42,5 °C (90,0 to 108,0 °F)
for surface measurements . . . . .	0 to 60 °C (32 to 140 °F)
Consumption . . . . .	≤ 50mV
Accuracy. . . . .	± 0,3 °C (0,54 °F)
Measuring distance . . . . .	5 cm - 15 cm (2,5 - 5,9 Zoll)
Automatic power-off . . . . .	after 7 seconds

### Accuracy of the non-contact infrared thermometer

32,0 to 35,9 °C / 93,2 to 96,6 °F . . . . . ± 0,3 °C (0,5 °F)

6,0 to 39,0 °C / 96,8 to 102,2 °F . . . . . ± 0,2 °C (0,4 °F)

39,0 to 42,5 °C / 102,2 to 108,5 °F . . . . . ± 0,3 °C (0,5 °F)

*Conforms to the ASTM standard E 1965-1998 (2003)*

The manufacturer reserves the right to change the technical specification of the model and to make modifications to the product as part of ongoing product development without prior notice.

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**TROTEC HealthCare GmbH • Grebbener Straße 7 • D-52525 Heinsberg**